

Next Week Glance

EXAM 3 INFO:

**Tues. March 28th
6:15 pm – 7:30 pm**

MON	TUE	WED	THU	FRI
3/27	3/28	3/29	3/30	3/31
Quiz 10	EXAM 3	PLQ 28 8 am	PLQ 29 R11 GQ	

Review INFO:

3/25 (Saturday) 1-3 pm, 101 Thomas, Practice Exam 3S

3/26 (Sunday) 1-2:30 pm, 304 Boucke, GSG

3/26 (Sunday) 3-5 pm, 102 Thomas, Practice Exam 3R

3/26 (Sunday) 6:30-8:30 pm, 105 Forum, TPE


3/27 (Monday) 3-5 pm, 117 Thomas, Practice Exam 3Q

Lesson 22-2 Weak Acid-Base Titrations

Lesson goals:

- Graphically describe the following titration curves; weak acid titrated with strong base, and weak base titrated with strong acid.
- If present, identify the equivalence point, half-equivalence point, and buffer region on a titration graph.
- Using tabulated K_a or K_b values and initial concentrations, calculate the pH at any point in any of the above titrations.
- Qualitatively predict the pH at the equivalence point of any of the above titrations. Describe the species present in solution in significant amounts at any point of the above titrations.
- Given information about the half-equivalence point, calculate the K_a or pK_a of the weak acid, or the K_b or pK_b of the weak base in a titration.
- Identify the titration curve of a polyprotic acid, and calculate the pH at any of the equivalence points.


Titration Curve for Weak Acid with a Strong Base




initial pH = 2.87;
higher than
strong acid

Equivalence Point;
pH > 7
only conjugate base
present


Half Equivalence Point;
[HX] = [X⁻]
pH = pK_a




1. Weak Acid ONLY
before OH⁻ addition
(pH; acid dissoc.)



2. BUFFER!
Half Equivalence Point
(pH; H-H)



3. Weak Base ONLY
Equivalence Point
(pH; base hydrolysis)




4. Excess Strong Base
Past Equivalence Point
(pH; based on xs OH⁻)

Lecture 29 Activity 1 – Weak Acid-Base Titrations

We will titrate 50 mL of 0.1 M CH_3COOH with varying amounts of 0.2 M NaOH .

1. Sketch the titration curve for this system on the axes to the right. At what pH range do you expect the equivalence point to be?
2. Write the **complete neutralization reaction** and the **net ionic** reaction under the graph. What are the spectator ions?
3. Calculate the pH before any NaOH is added. (For acetic acid, $K_a = 1.8 \times 10^{-5}$)



4. **Before doing titration calculations, ALWAYS CALCULATE THE MOLES OF REACTANT THAT WILL BE TITRATED.** How many total moles of H^+ will be titrated here when the endpoint is reached?

Lecture 29 Activity 3 – *pH at Equivalence Point*

We will titrate 50 mL of 0.1 M CH₃COOH with varying amounts of 0.2 M NaOH.

The equivalence point (or endpoint) is where:
Total Moles of base present = Total moles of H⁺ present

1. We will add 25 mL of 0.2 M NaOH solution to the same CH₃COOH solution. How many moles of base are being added?

	CH ₃ COOH	+ OH ⁻ →	CH ₃ COO ⁻	+ H ₂ O
I				

C

2. There are still the same number of initial moles of acid. Fill in the ICF table (Use moles or mmols!)

F				
---	--	--	--	--

3. What is the limiting reactant?

I				
---	--	--	--	--

4. What non-spectator species is left in solution after reaction? Classify it.

C

E				
---	--	--	--	--


5. What is the total volume?

6. What is the concentration of this species?

7. What is the pH of the solution at this point?
 ($K_a = 1.8 \times 10^{-5}$ for acetic acid.)







Acid Strength Affects Titration Curve Shape

For weaker acid,
higher initial pH,
flatter curve,
smaller ΔpH



Titration Curve of a Di- or Tri-protic Acid


Titrate 30 mL of 0.1 M H_3PO_4 solution with 0.1 M NaOH solution. $K_{a1} \gg K_{a2} \gg K_{a3}$



Titration of Bases Mirror Acid Titrations

Strong Base

Weak Base



Be able to predict **pH** at the **Equivalence Point**

- 1) Titrating a **strong acid** with a **strong base** (OR vice versa)

pH at equivalence = 7

- 2) Titrating a **weak acid** with a **strong base**

pH at equivalence > 7


★ *the weaker the acid, the smaller the pH change at equivalence*

- 3) Titrating a **weak base** with a **strong acid**

pH at equivalence < 7

The Indicator Must Change Color Within The pH Range at the Equivalence Point

methyl red (color change at pH 4.4 – 6.2)




Weak Acid:
indicator must
change color
between
pH = 7 to 10.8


Strong Acid:
indicator must
change color
between
pH = 2.8 to 10.8

Titration curves of NH_3 and NaOH are shown.
Which indicator is appropriate for the titration of NH_3 with HCl ?

methyl red (pH 4.4 – 6.2)



phenolphthalein (pH 8.0 – 9.8)



Indicator must change color between pH = 3 to 11 for strong base titration.
Indicator must change color between pH = 3 to 7 for weak base titration.

- A. methyl red only
- B. phenolphthalein only
- C. both methyl red and phenolphthalein are appropriate